

# EchoLUTION Tissue DNA Micro Kit – Protocol

for single-step purification of genomic DNA from tissue samples

This protocol has been developed for fresh, frozen and stabilized human and animal tissue samples including DNase-rich (e. g., spleen, liver, kidney) and lipid tissues (e.g. brain, fat).

## Materials and equipment needed

- 1 to 20 mg tissue per sample depending on tissue type.  
Sample input recommendation:

Tissue type	Generic	High DNA content (e.g., spleen, liver, kidney)	Low DNA content (e.g., muscle, cartilage)
Amount	10 mg	5 mg	20 mg

- Microcentrifuge with rotor for 1.5 and 2 ml reaction tubes.  
**Important:** Switch centrifuge to *relative centrifugal force, rcf* ( $\times g^*$ ); if this is not possible, please use formula below\* to calculate the conversion of rounds per minute (rpm) into rcf.
- Thermal shaker with agitation (for fastest performance), capable of heating to 60°C and 80°C. Alternatively: Heating Block or heat chamber
- Vortexer
- Pipets for 10  $\mu$ l and 200  $\mu$ l scales, corresponding pipet tips
- One reaction tube (1.5 ml) per sample for the lysis step (preferably safe-lock).
- One reaction tube (2 ml) per sample for column preparation.
- One reaction tube (1.5 ml) per sample for elution and collection of the purified DNA.
- Ceramic Blade Scalpels. BioEcho product no. 050-002-001
- For fastest procedure (PROTOCOL 1): Cap Puncher (BioEcho product no. 050-001-001)

## Preparation before starting

- Heat the thermal shaker or thermo block to 60°C
- Set the microcentrifuge to 1,000  $\times g^*$
- Important:** Switch to *relative centrifugal force, rcf* ( $\times g^*$ , not rpm)

## PROTOCOL 1: Purification using the Cap Puncher



### Lysis

- Transfer tissue sample to the bottom of a 1.5 ml reaction tube (preferably safe-lock) while the tube is cooled on ice (or cooling block) to avoid DNA degradation during sample loading.

#### Note:

- If possible, cut tissue into small pieces to speed up lysis.
- For stabilized tissue samples briefly rinse with water to remove traces of stabilization solution before adding samples to the reaction tube.

- Add 90  $\mu$ l Tissue Lysis Buffer and 5  $\mu$ l TurboLyse T Protease to each tissue sample. Mix by flicking or vortexing.

If working with more than two samples, prepare a Lysis Master Mix with 10% excess volume for the number of tissue samples (see table).

### Lysis Master Mix:

No of samples	1	6 (+10%)	12 (+10%)	Yours
<b>P</b> Tissue Lysis Buffer ( $\mu$ l)	90	600	1,190	
<b>LB</b> TurboLyse T Protease ( $\mu$ l)	5	35	70	
Final volume ( $\mu$ l)	95	635	1,260	

Add 95  $\mu$ l of the Lysis Master Mix to each tissue sample.

- Place the reaction tube(s) in the thermal shaker and incubate at 60°C for 30 min with max. agitation (for 60 min if agitation is not feasible, in this case pulse-vortex 3 times during lysis).

**Note:** If samples are not completely lysed after the time period described above, continue with the next step. Residual cellular debris will not interfere with the purification performance.

**Note:** For some tissue types, lysis is already complete after 15 min. Step 3 may be shortened accordingly.

Meanwhile during lysis, proceed with step 5, "Column preparation".

- After incubation at 60°C, increase the temperature to 80°C and incubate for additional 10 min with max. agitation.

### Column preparation (during step 3. and 4.)


- Vortex the EchoLUTION Spin Column briefly and place into a 2 ml reaction tube. Let stand for 10-20 min.
- Use of the Cap Puncher (scan QR code to watch a video): Punch a hole into the column cap and lift the column together with the Cap Puncher out of the 2 ml collection tube. Snap off bottom closure of the column and detach the Cap Puncher by twisting while pulling out. Place the punched spin column back into the 2 ml reaction tube.
- Centrifuge for 1 min at 1,000  $\times g^*$ . Discard the flow-through volume ("void volume") collected in the 2 ml reaction tube.
- Place the prepared EchoLUTION Spin Column into a new 1.5 ml reaction tube for elution of the purified DNA and place back into the rack. Continue with "Purification of DNA".

### Purification of DNA

- Optional:** After having performed step 4, add 1  $\mu$ l RNase **R** to each lysed sample and vortex 3 sec. Incubate for 2 min at room temperature to remove RNA traces.
- After having performed step 4, add 10  $\mu$ l Clearing Solution T **CS** each sample. Vortex 3 sec. The sample becomes cloudy.
- Centrifuge for 2 min at full speed.
- Transfer lysis supernatant (max. 100  $\mu$ l) containing the DNA onto the prepared EchoLUTION Spin Column from step 8 as illustrated: Insert pipet tip vertically through the hole in the column cap and pipet the sample slowly (~5 sec) into the column.

#### Note:

- During loading of lysate, do not touch the resin bed!
  - Residual sample particles may be loaded and will not interfere with purification.
- Centrifuge for 1 min at 1000  $\times g^*$ . The purified DNA elutes into the 1.5 ml elution tube and can be immediately applied in downstream applications.



## PROTOCOL 2: Purification without a Cap Puncher

### Lysis


1. Perform steps 1-4 from PROTOCOL 1.

### Column preparation

5. Vortex the EchoLUTION Spin Column briefly and place into a 2 ml reaction tube. Let stand for 10-20 min (recommended to be done during lysis step).
6. Loosen the screw cap of the spin column half a turn and snap off the bottom closure. **Important:** Do not re-close the screw cap of the spin column. The screw cap must stay loosened half a turn to avoid generation of a vacuum. Place the column back into the 2 ml collection tube and both into the centrifuge.
7. Centrifuge for 1 min at 1,000 x g\*. Discard the 2 ml reaction tube containing the column buffer.
8. Place the prepared spin column into a new 1.5 ml reaction tube for elution of the sample DNA and place back into the rack. Continue with "Purification of DNA".

### Purification of DNA

9. After having performed step 4, add 1 µl RNase **(R)** each lysed sample and vortex for 3 sec. Let stand for 2 min at room temperature to remove RNA traces.
10. Add 10 µl Clearing Solution T **(CS)** each sample. Vortex 3 sec. The sample will become cloudy.
11. Centrifuge for 2 min at maximum speed.
12. Transfer lysate supernatant (max. 100 µl) containing the DNA onto the prepared column from step 8 as illustrated:



Open cap and pipet the sample slowly (~5 sec) onto the center of the resin bed of the prepared spin column. Close screw cap and loosen again half a turn.

**Important:** Do not re-close the screw cap of the spin column completely!

### Note:

- During loading of lysate, do not touch the resin bed with the pipet tip!
  - Residual sample particles may be loaded and will not interfere with purification.
13. Centrifuge 1 min at 1,000 x g\*. The purified DNA elutes into the 1.5 ml elution tube and can be immediately applied in downstream applications.

### Product use limitation

The EchoLUTION Tissue DNA Micro Kit is for research use only. It is not registered or authorized to be used for diagnosis, prevention or treatment of a disease.

\* Most centrifuges offer the choice between rpm and g-force (rcf); if not, calculate the rpm matching the g-force using the formula:  $\text{rpm} = 1,000 \times \sqrt{\frac{g}{1.12 \times r}}$ , where r = radius of rotor in mm. and g the required g-force. E. g., with a radius of 150 mm, the corresponding rpm to 1,000 x g is approx. 2,400 rpm.

# EchoLUTION Tissue DNA Micro Kit

for single-step purification of genomic DNA from tissue samples

Product no. (rxn's)	010-002-010 (10)	010-002-050 (50)	010-002-250 (250)
Kit contents	Tissue Lysis Buffer, TurboLyse T Protease, RNase A Tissue, Clearing Solution T, 1x Tris Buffer, EchoLUTION Spin Columns		
Related products	Cap Puncher 050-001-001	Ceramic Blade Scalpels 050-002-00x	

## Quick PROTOCOL (please read protocol first)

### Lysis


- Transfer 10 mg tissue to reaction tube
- Add 90 µl **(LB)**
- Add 5 µl **(P)**, vortex briefly
- Incubate 30 min/60°C, then 10 min/80°C, max. agitation
- Add 1 µl **(R)**, vortex shortly and incubate 2 min at room temperature
- Add 10 µl **(CS)** and vortex shortly
- Centrifuge 2 min at max. speed

### Column preparation (during 60°C and 80°C incubation)

- Vortex EchoLUTION spin column and place in a 2 ml tube. Let stand for 10 min
- Punch a hole in the cap with the Cap Puncher, and break off bottom closure (scan QR code to watch a video)
- Place spin column back into 2 ml tube
- Centrifuge 1 min at 1,000 x g\* to elute column buffer
- Place column in a 1.5 ml tube

### Purification of DNA

- Transfer lysate supernatant (max. 100 µl) by pipetting slowly through cap hole (scan QR code to watch a video)
- Centrifuge 1 min at 1,000 x g\* to elute DNA into Elution tube
- Eluted DNA is ready to use



2020/01/23



BioEcho Life Sciences GmbH  
Nattermannallee 1  
50829 Köln (Cologne)/Germany

Phone: +49 (0) 221-99 88 97-0  
E-Mail: [contact@bioecho.de](mailto:contact@bioecho.de)  
[www.bioecho.de](http://www.bioecho.de)