

# EchoCLEAN DNA CleanUp Kit – Protocol

for single-step depletion of impurities and DNA fractions < 50bp from DNA solutions

This protocol has been developed to deplete impurities (e. g., salts, peptides, nucleotides and DNA fragments and primers < 50 bp) from DNA solutions.

## Materials and equipment needed

- 80 to 110 µl of DNA sample. For DNA samples less than 80 µl, dilute with Tris buffer to a minimum of 80 µl
- Microcentrifuge with rotor for 1.5 and 2 ml reaction tubes  
**Important:** Switch centrifuge to *relative centrifugal force, rcf* ( $x g^*$ ); if this is not possible please use formula below\* to calculate the conversion of round per minute (rpm) into rcf.
- Vortexer
- One reaction tube (2 ml) per sample for column preparation
- One reaction tube (1.5 ml) per sample for elution and collection of the purified DNA
- Pipets for 10 µl and 200 µl scale, corresponding pipet tips
- For fastest procedure (PROTOCOL 1): Cap Puncher (BioEcho product no. 050-001-001)

## Preparation before starting

- Set the microcentrifuge to  $1000 x g^*$
- **Important:** Switch to *relative centrifugal force, rcf* ( $x g^*$ , not rpm)

## Product use limitation

The EchoCLEAN DNA CleanUp Kits are for research use only. They have not been registered or authorized to be used for diagnosis, prevention or treatment of a disease.

**Note:** The depletion of nucleic acids of a specific length is a function of loaded concentration and amount of fragment to be depleted. In certain cases, depletion may not be quantitative

\* Most centrifuges offer the choice between rpm and g-force (rcf); if not, calculate the rpm matching the g-force using the formula:  $rpm = 1,000 \times \sqrt{\frac{g}{1.12 \times r}}$ , where r = radius of rotor in mm. and g the required g-force.  
E. g., with a radius of 150 mm, the corresponding rpm to  $1,000 x g$  is approx. 2,400 rpm.

## PROTOCOL 1: Purification using the Cap Puncher




BioEcho Cap Puncher

### Column preparation

1. Vortex the EchoCLEAN DNA Spin Column briefly and place into a 2 ml reaction tube. Let stand for about 5 min.
2. Use of the Cap Puncher (scan QR code to watch a video): Punch a hole into the column cap and lift the column together with the Cap Puncher out of the 2 ml collection tube. Snap off bottom closure of the column and detach the Cap Puncher by twisting while pulling out.  
Place the punched spin column back into the 2 ml reaction tube
3. Centrifuge for 1 min at  $1000 x g^*$ . Discard the flow-through volume (“void volume”) collected in the 2 ml reaction tube.
4. Place the prepared spin column into a new 1.5 ml reaction tube for elution of the purified DNA and place back into the rack.  
Continue with “Clean up of DNA”.

### Clean up of DNA


5. Transfer 80 – 110 µl of the DNA sample to the prepared EchoCLEAN Spin Column from step 4 as illustrated below:



Insert pipet tip vertically through the hole in the column cap and pipet the sample slowly (~5 sec) into the column.

### Note:

- Do not punch pipette tip into the resin bed during loading of lysate!
6. Centrifuge 1 min at  $1000 x g^*$ . The purified genomic DNA elutes into the 1.5 ml elution tube and can be immediately applied in downstream applications.

The eluted DNA can be used immediately or stored at  $4^{\circ}\text{C}$  or for long-term storage at  $-20^{\circ}\text{C}$ . For spectrophotometric analysis, use 1x Tris Buffer  supplied with the kit as blank.

# EchoCLEAN DNA CleanUp Kit

for single-step depletion of impurities and DNA fractions < 50bp from  
DNA solutions


## PROTOCOL 2: Purification without a Cap Puncher

### Column preparation

1. Vortex the EchoCLEAN DNA Spin Column briefly and place into a 2 ml reaction tube. Let stand for 5 min.
2. Loosen the screw cap of the spin column half a turn and snap off the bottom closure.  
**Important:** Do not close the screw cap of the spin column. The screw cap must stay loosened half a turn to avoid generation of a vacuum. Place the column back into the 2 ml collection tube and both into the centrifuge.
3. Centrifuge for 1 min at  $1000 \times g^*$ . Discard the flow-through volume (“void volume”) collected in the 2 ml reaction tube.
4. Place the prepared spin column into a new 1.5 ml reaction tube for elution of the sample DNA and place back into the rack.  
Continue with “Clean up of DNA”.

### Clean up of DNA

5. Transfer 80 – 110  $\mu$ l of the DNA sample to the prepared EchoCLEAN Spin Column from step 4 as illustrated below




Open cap and pipet the sample slowly (~5 sec) onto the center of the resin bed of the prepared spin column. Close screw cap and loosen again half a turn.

**Important:** Do not close the screw cap of the spin column tightly!

#### Note:

- Do not punch pipette tip into the resin bed during loading of lysate!
6. Centrifuge 1 min at  $1000 \times g^*$ . The purified genomic DNA elutes into the 1.5 ml elution tube and can be immediately applied in downstream applications

The eluted DNA can be used immediately or stored at  $4^\circ\text{C}$  or for long-term storage at  $-20^\circ\text{C}$ . For spectrophotometric analysis, use 1x Tris Buffer  supplied with the kit as blank.

Product no. (rxn's)	020-002-030-010 (10)	020-002-030-050 (50)	020-002-030-250 (250)
Kit contents	EchoCLEAN DNA CleanUp Spin Columns, 1x Tris Buffer		


### Quick PROTOCOL (please read protocol first)

#### Column preparation

- Vortex EchoCLEAN spin column and place in a 2 ml tube  
Let stand for 10 min
- Punch a hole in the cap with the Cap Puncher, and break off bottom closure (scan QR code to watch a video)
- Place spin column back into 2 ml tube
- Centrifuge 1 min at  $1000 \times g^*$  to elute column buffer
- Place column in a 1.5 ml tube

#### Clean up of DNA

- Transfer DNA sample (max. 110  $\mu$ l) by pipetting slowly through cap hole (scan QR code to watch a video)
- Centrifuge 1 min at  $1,000 \times g^*$  to elute DNA into Elution tube
- Eluted DNA is ready to use



2020/01/23



BioECHO Life Sciences GmbH  
Nattermannallee 1  
50829 Köln (Cologne)/Germany

Phone: +49 (0) 221-99 88 97-0  
E-Mail: [contact@bioecho.de](mailto:contact@bioecho.de)  
[www.bioecho.de](http://www.bioecho.de)