


Selection Guide EchoCLEAN CleanUp Kits

By Impurity

Kit \ Impurity	Organic solvents				DNA		Salts				Dyes		Other ¹
	Phenol	Trizol	Chloroform	Ethanol	Primer	dNTPs	Chao-trophs ²	Salts	SDS	Na-Azide	Indigo	Gel loading	Precipitates
EchoCLEAN DNA CleanUp Kit (for DNA > 50 bp)				✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
EchoCLEAN Organic Solvent DNA CleanUp Kit	✓	✓	✓	✓				✓					✓
EchoCLEAN RNA CleanUp Kit ³	✓	✓	✓	✓			✓	✓					✓

New!
New!



By Application

Kit \ Application	Phenol-based Extraction ⁴	Desalting of Nucleic Acids	Post Bind-Wash-Elute Extraction	PCR CleanUp	Enzymatic Reaction CleanUp	Buffer Exchange ⁵
EchoCLEAN DNA CleanUp Kit (for DNA >50 bp)		✓	✓	✓	✓	✓
EchoCLEAN Organic Solvent DNA CleanUp Kit	✓	✓				✓
EchoCLEAN RNA CleanUp Kit ³	✓	✓	✓			✓

¹ CleanUp of humic acids, please inquire
² Gu-HCl; GTC
³ does not recover miRNAs, tRNAs, small RNA molecules
⁴ Trizol; Phenol/Chloroform; based – extractions
⁵ Buffering to low mM Tris, pH 8.2

convenience • performance • sustainability